

Sterile Media Bottles


PET Bottle with HDPE closure

- **☑** Radiation Sterile
- **✓** Non-Cytotoxic **☑** Non-Pyrogenic
- **☑** DNase / RNase free
- **☑** Conforming to USP ClassVI
- ☑ Animal derived component free



Sterile Media Bottles

- > An inexpensive alternative to glass bottles for sampling, storage and shipping of culture media, API and sera with a guaranteed leak-prood seal and exellent gas barrier proper-
- Light weight, shatter resistant and space saving square
- Transparent bottles made of non-cytotoxic, non-mutagenic, non-pyrogenic and Animal Derived Component free virgin resins to ensure content purity.
- > Keeps Products safe by meeting the latest biological testing standards including USP Class VI, USP37<661>, and ISO
- > Eliminate costly washing, depyrogenation and autoclaving



Square PET Media Bottles, with HDPE closure					
Cat. No.	Capacity	Closure Size	Quantity per tray/case		
581610	125ml	38-430	24/96		
581611	250ml	38-430	30/60		
581612	650ml	38-430	24/48		
581613	1000ml	38-430	12/24		

Square PETG Media Bottles, with White HDPE closure					
Cat. No.	Capacity	Closure Size	Quantity per tray/case		
581410	30ml	20-415	24/96		
581411	60ml	24-415	24/96		
581412	125ml	38-430	24/48		
581413	250ml	38-430	24/48		
581414	500ml	38-430	12/24		
581415	1000ml	38-430	12/24		
581416	2000ml	53B	6/12		

^{**}Square PETG Bottle Design - Patent Pending

Round PETG Diagnostic Bottles, with White HDPE closure					
Cat. No.	Capacity	Closure Size	Quantity per tray/case		
581701	5ml	20-415	20/100		
581702	10ml	20-415	20/100		
581703	20ml	20-415	20/100		

Special Feature:

- Square Bottles are graduated
- Square PETG Bottles features heat-shrink band around closure and neck.
- Maximum Use Temp for PET bottles is 65°C
- Maximum use Temp for PETG bottle is 70°C



TARSONS PRODUCTS LIMITED

902, Martin Burn Business Park, BP-3, Saltlake, Sector - V, Kolkata - 700091 Tel: +91 033 - 35220300, e-mail: info@tarsons.com Website: www.tarsons.com



TarsonsProductsPvtLtd in tarsons-products-pvt-ltd otarsons.products TarsonsProduct





