

Sterile Sampling Tools

Sampling Spoons, Scoops, Spatulas, and Forceps

- ▶ Radiation Sterile
- ▶ Conforming to USP Class VI, FDA 21 CFR
- ▶ Disposable
- ▶ Non-Cytotoxic
- ▶ Individually packed



Special features:

- ▶ Long handles make them a perfect option for easy collection and aseptic sampling of your precious samples from hard-to-reach areas
- ▶ Safe for food, drug and cosmetic use
- ▶ Resistant to a temperature range of -40°C to 80°C for HIPS & 0°C to 80°C for PS
- ▶ V-shaped spatulas are ideal for piercing plastic-lined paper bags as well as sampling of granulated/semi solid material
- ▶ Corrugated handles of forceps provide a secure grip
- ▶ Reduces cross contamination

Sterile Sampling Tools

Ordering information:

Sterile Individually Packed Long Handle Sampling Spoon, HIPS

Code	Capacity	Handle Length	Packing
922121	1.25 ml	155.4 mm (6")	200
922122	2.5 ml	149.5 mm (5.8")	200
922123	5 ml	143.5 mm (5.6")	200
922124	15 ml	129.5 mm (5.1")	200
922111	1.25 ml	155.4 mm (6")	10
922112	2.5 ml	149.5 mm (5.8")	10
922113	5 ml	143.5 mm (5.6")	10
922114	15 ml	129.5 mm (5.1")	10

Note: Overall length of all spoons is 178mm (7")

Sterile Individually Packed Long Handle Sampling Scoop, HIPS

Code	Capacity	Handle Length	Overall Length	Packing
923111	60 ml	126.5 mm (5")	204 mm (8")	100
923112	125 ml	126.5 mm (5")	220 mm (8.6")	100
923113	250 ml	88 mm (3.5")	214.3 mm (8.4")	100
923101	60 ml	126.5 mm (5")	204 mm (8")	10
923102	125 ml	126.5 mm (5")	220 mm (8.6")	10
923103	250 ml	88 mm (3.5")	214.3 mm (8.4")	10

Note: Overall length of all spatulas is approximately 229mm (9")

Sterile Individually Packed V-shaped Sampling Spatulas, HIPS

Code	Capacity	Blade Length	Handle Length	Packing
921130	15 ml	152.4 mm (6")	76.2 mm (3")	100
921120	15 ml	152.4 mm (6")	76.2 mm (3")	10

Sterile Individually Packed Disposable Forceps, PS

Code	Size	Length	Packing
486001	124.2 X 21.0 X 8.0 mm	4.8"	100



TARSONS PRODUCTS LIMITED

902, Martin Burn Business Park,
BP-3, Saltlake, Sector - V, Kolkata - 700091
Tel: +91 033 - 35220300, e-mail: info@tarsons.com
Website: www.tarsons.com

