


# SCREW CAP MICRO TUBES

*~Ideal for transport and storage of enzymes, oligo-nucleotides, buffers etc.*



## Ensure Reliable Sealing Security



For more information, visit us at: [www.tarsons.com](http://www.tarsons.com)

# SCREW CAP MICRO TUBES



Cap with EPDM 'O' ring



Tube with Knurl



Tube with Looped Cap

## Screw Cap Tubes - Non-Sterile (without Cap)

Code	Capacity mL	Description	Packing unit/case
530010	0.5	Self standing with knurl	500/5000
530011	0.5	Self standing without knurl	500/5000
530020	0.5	Conical with knurl	500/5000
530030	1.5	Self standing with knurl	500/5000
530040	1.5	Conical with knurl	500/5000
530041	1.5	Conical without knurl	500/5000
530050	2.0	Self standing with knurl	500/5000
530051	2.0	Self standing without knurl	500/5000
530060	2.0	Conical with knurl	500/5000

## Screw Cap Tubes - Sterile (without Cap)

Code	Capacity mL	Description	Packing unit/case
S530010	0.5	Self standing with knurl	500/5000
S530011	0.5	Self standing without knurl	500/5000
S530020	0.5	Conical with knurl	500/5000
S530030	1.5	Self standing with knurl	500/5000
S530040	1.5	Conical with knurl	500/5000
S530041	1.5	Conical without knurl	500/5000
S530050	2.0	Self standing with knurl	500/5000
S530051	2.0	Self standing without knurl	500/5000
S530060	2.0	Conical with knurl	500/5000

## Screw Cap Tubes - Amber (without Cap)

Code	Capacity mL	Description	Packing unit/case
A530010	0.5	Self standing with knurl	500/5000
A530011	0.5	Self standing without knurl	500/5000
A530020	0.5	Conical with knurl	500/5000
A530030	1.5	Self standing with knurl	500/5000
A530040	1.5	Conical with knurl	500/5000
A530041	1.5	Conical without knurl	500/5000
A530050	2.0	Self standing with knurl	500/5000
A530051	2.0	Self standing without knurl	500/5000
A530060	2.0	Conical with knurl	500/5000

## CAP for Screw Cap Tubes (Sterile/Non-Sterile)

Code	Description	Packing unit/case
530000	Cap with EPDM 'O' ring Non-Sterile	500/5000
S530000	Cap with EPDM 'O' ring Sterile	500/5000

## Screw Cap Tubes with Looped Cap unassembled Non-Sterile

Code	Capacity mL	Description	Packing unit/case
532010	0.5	Self standing with knurl	500/5000
532030	1.5	Self standing with knurl	500/5000
532040	1.5	Conical with knurl	500/5000
532050	2.0	Self standing with knurl	500/5000
532060	2.0	Conical with knurl	500/5000

## Screw Cap Tubes with Looped Cap Sterile

Code	Capacity mL	Description	Packing unit/case
S532010	0.5	Self standing with knurl	500/5000
S532030	1.5	Self standing with knurl	500/5000
S532040	1.5	Conical with knurl	500/5000
S532050	2.0	Self standing with knurl	500/5000
S532060	2.0	Conical with knurl	500/5000

Add suffix N= Natural, W= White, R= Red, B= Blue, Y= Yellow, G= Green, Z= Amber for color cap.



Conical with knurl - 0.5 mL



Conical without knurl-1.5 mL



Self standing without knurl- 2 mL



Tube with unassembled Looped Cap (Non-Sterile)



Tube with Looped Cap (Sterile)- 2 mL



Self standing with knurl - 2 mL

- Made of high grade Polypropylene material, conforming to USP class VI.
- Free of detectable DNase, RNase and Pyrogen / Endotoxin.
- Free of detectable Human DNA and PCR inhibitors.
- Heavy metal free EPDM 'O' ring.
- Leak-proof screw cap with O-ring.
- Max RCF - 20,000 X g.
- Manufactured in ISO Class-8 Cleanroom.

Note: Not recommended for liquid Nitrogen submersion storage. Short term vapour phase and Quick Freeze acceptable.



TarsonsProductsLimited



tarsons-products-ltd



TarsonsProduct



tarsons.products