

Sterile Square Biocanisters with white HDPE Closure



Available in
5L, 10L & 20L



- ✓ Leakproof & Graduated
- ✓ Non-cytotoxic & Free of detectable pyrogens
- ✓ Available in PC, PET & PETG
- ✓ 3-layer packaging for cleanroom use

New Launch

Ideal choice for safe & cost-effective storage



Sterile Square Biocanisters

~ Ordering information:

Sterile Square PC Biocanister, with white HDPE Closure

Material conforming to USP Class VI. Non cytotoxic and free of detectable pyrogens. Graduated. Individually packed. Ideal for safe storage and transport of vaccines, biopharmaceuticals, culture media, and other biotech materials in life science industries. Withstands -80°C to 121°C. Moulded in graduations with 500 mL/0.5-gallon increments in 5L; 500mL/1-gallon increments in 10L and 1 L/1-gallon increments in 20L.

Code	Capacity (L)	Closure size (mm)	Packing
581407	5	63	6
581408	10	63	4
581409	20	63	2

Sterile Square PET Biocanister, with white HDPE Closure

Material conforming to USP Class VI. Non cytotoxic and free of detectable pyrogens. Graduated. Individually packed. Ideal for safe storage and transport of vaccines, biopharmaceuticals, culture media, and other biotech materials in life science industries. Withstands -40°C to 65°C. Moulded in graduations with 500 mL/1-gallon increments in 10L and 1 L/ 1-gallon increments in 20L.

Code	Capacity (L)	Closure size (mm)	Packing
581418	10	63	4
581419	20	63	2

Sterile Square PETG Biocanister, with white HDPE Closure

Material conforming to USP Class VI. Non cytotoxic and free of detectable pyrogens. Graduated. Individually packed. Ideal for safe storage and transport of vaccines, biopharmaceuticals, culture media, and other biotech materials in life science industries. Withstands -40°C to 70°C. Moulded in graduations with 500-mL and 0.5-gallon.

Code	Capacity (L)	Closure size (mm)	Packing
581417	5	63	6

Note: Gamma stable - Non sterile biocanisters are also available.



~Technical specifications may be changed without prior notice.



TARSONS PRODUCTS LIMITED

902, Martin Burn Business Park,
BP-3, Saltlake, Sector - V, Kolkata - 700091
Tel: +91 033 - 35220300, e-mail: info@tarsons.com
Website: www.tarsons.com

