


Petri Dish



CONTACT PLATE & ACCESSORIES



- ✓ High Optical clarity
- ✓ Ultra flat base
- ✓ Easy to stack
- ✓ Ideal for use in automatic plate pourers
- ✓ Aseptically manufactured in ISO Class 8 Clean Room



Culture micro-organisms with confidence




Disposable Radiation Sterile Petri Dishes Triple Vent, PS

Material : PS conforming to US FDA 21 CFR

STERILE R

Vented dishes have elevated lid which allows for good gas exchange and encourage evaporation. This design is useful when the poured dishes are to be used as soon as possible. The design of the dish maintains sterility because the particles would have to go up and over the dish wall to get inside which is very rare in normal airflow.

Code	Size (dia x h) mm	Packing
460035	35 x 15	500 pcs individual packing
460038	35 x 15	500 pcs packed in sleeve of 10
460062	60 x 15	250 pcs individual packing
460061	60 x 15	360 pcs packed in sleeve of 15
460090	90 x 14	450 pcs individual packing
460091	90 x 14	480 pcs packed in sleeve of 10
460095	90 x 14	480 pcs packed in sleeve of 20
460140	140 x 16	136 pcs individual packing



Disposable Radiation Sterile Petri Dishes Non Vented

Material : PS conforming to US FDA 21 CFR

STERILE R

Non vented dishes have lid which fits flatly on the base resulting in less potential for external contamination and significantly reduced evaporation rate.

Code	Size (dia x h) mm	Packing
460020	90 x 14	450 pcs individual packing
460030	90 x 14	480 pcs packed in sleeve of 10
460040	90 x 14	480 pcs packed in sleeve of 20

Disposable Aseptic Petridish, PS

Material : PS conforming to US FDA 21 CFR

These Petri Dishes are manufactured in ISO class 8 clean room under aseptic conditions which ensures no contact with any contaminants. Aseptically produced Petridishes are particulate free and have the added advantage of being free of chemical residues produced by ETO which causes false negative reading. It also eliminates discolouration and stress caused by Radiation Sterilization.

Code	Size (dia x h) mm	Packing
Triple Vent		
460051	90 x 14	450 pcs individual packed
460092	90 x 14	500 pcs packed in sleeve of 20
Non Vented		
460050	90 x 14	500 pcs packed in sleeve of 20



Petri Dish Autoclavable, PP

Material : PP conforming to US FDA 21 CFR, Autoclavable

Code	Size (dia x h) mm	Packing
462010	50 x 14	12 pcs
462020	75 x 14	12 pcs
462030	100 x 14	12 pcs
462040	125 x 20	12 pcs

Petri Dish Autoclavable, PMP

Material : PMP (TPX) conforming to US FDA 21 CFR, Autoclavable

Code	Size (dia x h) mm	Packing
461010	50 x 14	12 pcs
461020	75 x 14	12 pcs
461030	100 x 14	12 pcs
461040	125 x 20	12 pcs





Disposable Radiation Sterile Contact Plate, PS

Material : PS conforming to US FDA 21 CFR

STERILE R

Contact plates are used in routine hygiene monitoring of surfaces. Vent aids gas exchange and numbered grid on the base facilitates colony counting.

Code	Size (dia x h) mm	Base Type	Packing
460065	55 x 14	Griplock Flat	400 pcs packed in sleeve of 10
460067	55 x 14	Griplock Dome	360 pcs packed in sleeve of 10
460070	55 x 14	Dome	360 pcs packed in sleeve of 10




Disposable Radiation Sterile Petri Dish Square 125 mm x 125 mm, PS

Material : Clear PS conforming to US FDA 21 CFR.

STERILE R

Stackable and optimized for use on automated platform. Used in colony picking and antibiogram applications.

Code	Size (l x h)mm	Type	Packing
460120	120 X 17	Vented	240 pcs packed in sleeve of 10
460121	120 X 17	Non-vented	



Disposable Radiation Sterile Petri Dish Triple Vent Extra Depth, PS

Material : PS conforming to US FDA 21 CFR

STERILE R

The extra depth is an advantage for plant tissue culture experiment.

Code	(dia x h) mm	Packing
460094	90 x 20	360 pcs individual packing
460096	90 x 20	360 pcs packed in sleeve of 15




Disposable Sterile Soft Loop , PS

STERILE | R

This non-toxic disposable loop is accurate to semi-quantitative standard for handling sample. The surface is treated to enhance adhesion of droplets.

Code	Description	Packing
920051	Soft Loop - 1 μ l	100 pcs in pouch of 10
920061	Soft Loop - 10 μ l	
920071	Inoculation Needle	100 pcs in pouch of 20



Disposable Sterile L Shaped Spreader , PS

STERILE | R

This sterile, smooth spreader enables even spreading of liquid sample across the surface of agar plate without gouging or cutting the medium.

Code	Description	Packing
920082	L shaped spreader	100 pcs individual packing



Rack for Petri Dish, PMMA/PC

Code	Places	Type	Packing
270090	60	For 90 mm Petri Dish	2 pcs
270060	56	For 60 mm Petri Dish	2 pcs
270070	56	For Contact Plate	2 pcs



Incubation Tray, RPP

Code	Packing
525190	4 pcs

~Technical specifications may be changed without prior notice



TARSONS PRODUCTS LIMITED

902, Martin Burn Business Park,
BP-3, Saltlake, Sector - V, Kolkata - 700091
Tel: +91 033 - 35220300, e-mail: info@tarsons.com
Website: www.tarsons.com