

Sterile Sampling Tools



Sampling Spoon, Scoop, Spatula, and Forcep

- ✓ Sterile / Autoclavable
- ✓ Conforming to USP Class VI, FDA 21 CFR
- ✓ Disposable
- ✓ Non-Cytotoxic
- ✓ Individually packed



Special features:

- Perfect option for easy collection and aseptic sampling of precious samples
- Safe for food, drug and cosmetic use
- Resistant to a temperature range of -40°C to 80°C for HIPS, 0°C to 80°C for PS & 0°C to 121°C for PP
- V-shaped spatulas are ideal for piercing plastic-lined paper bags as well as sampling of granulated/semi solid material
- Forceps come with Serrated & Non-Serrated Blunt Tips
- Scissor-type forceps feature a scissor type handle with double action hinges that provide extra leverage & long life
- Scissor-type forceps are Non-Sterile, Autoclavable & Reusable with serration to grip objects tightly

Ideal choice for safe & aseptic sampling



Sterile Sampling Tools

Ordering information:

Sterile Individually Packed Long Handle Sampling Spoon, HIPS

Code	Capacity	Handle Length	Packing
922121	1.25 ml	155.4 mm (6")	200
922122	2.5 ml	149.5 mm (5.8")	200
922123	5 ml	143.5 mm (5.6")	200
922124	15 ml	129.5 mm (5.1")	200
922111	1.25 ml	155.4 mm (6")	10
922112	2.5 ml	149.5 mm (5.8")	10
922113	5 ml	143.5 mm (5.6")	10
922114	15 ml	129.5 mm (5.1")	10

Note: Overall length of all spoons is 178mm (7")

Sterile Individually Packed Long Handle Sampling Scoop, HIPS

Code	Capacity	Handle Length	Overall Length	Packing
923111	60 ml	126.5 mm (5")	204 mm (8")	100
923112	125 ml	126.5 mm (5")	220 mm (8.6")	100
923113	250 ml	88 mm (3.5")	214.3 mm (8.4")	100
923101	60 ml	126.5 mm (5")	204 mm (8")	10
923102	125 ml	126.5 mm (5")	220 mm (8.6")	10
923103	250 ml	88 mm (3.5")	214.3 mm (8.4")	10

Sterile Individually Packed V-shaped Sampling Spatula, HIPS

Code	Capacity	Blade Length	Handle Length	Packing
921130	15 ml	152.4 mm (6")	76.2 mm (3")	100
921120	15 ml	152.4 mm (6")	76.2 mm (3")	10

Note: Overall length of all spatulas is approximately 229mm (9")

Sterile Individually Packed Disposable Tweezer-type Forcep, PS

Code	Type	Length	Packing
486001	Non-Serrated	124.2 mm (4.88")	100
486002	Serrated	124.2 mm (4.88")	100

Scissor-type Forcep, PP

Code	Type	Length	Packing
486000	Serrated	122 mm (4.8")	12

Note: Non-Sterile items are also available

~ Technical specifications may be changed without prior notice



TARSONS PRODUCTS LIMITED

902, Martin Burn Business Park,
BP-3, Saltlake, Sector - V, Kolkata - 700091
Tel: +91 033 - 35220300, e-mail: info@tarsons.com
Website: www.tarsons.com