

Roller Bottles



~Roll into Superior Cell Growth with our High performance ROLLER BOTTLES



- Made from the highest quality, optically clear, virgin Polystyrene (PS) and Polyethylene Terephthalate (PET)
- Comes with virgin High-density Polyethylene (HDPE) caps
- Available in both Smooth & Ribbed surface formats
- Cell culture treated, for consistent, reliable cell attachment
- Free of detectable DNase, RNase, and human DNA
- Endotoxin testing conducted as per USP 85 with tolerance level of <0.01 EU/mL
- Sterilised by gamma irradiation, validated according to ISO 11137, which meets SAL 10⁻⁶
- All bottles are graduated with printed LOT number & Expiry Date
- Certified non-cytotoxic and non-pyrogenic
- Manufactured in ISO 13485 certified Class 7 Cleanrooms

For Consistent & Reliable Cell attachment



Sterile Roller Bottles, Smooth Wall, PS

Made from optically clear, virgin, [Polystyrene](#), these bottles are available in a [smooth surface](#) with a standard growth area of [850 cm²](#) and [1700 cm² \(2X\)](#). Available either with seal screw caps or vent screw caps having 0.2µ PTFE membrane closures that ensure contamination-free cultivation.



Code	Surface Area	Nominal Capacity	Dimensions (od x h)	Cap Style	Cap Colour	Surface	Pack/Case
447021	850 cm ²	2000 mL	122 x 272 mm	Seal Cap	Blue	TC Treated	1/12
447031	850 cm ²	2000 mL	122 x 272 mm	Vent Cap	Blue	TC Treated	1/12
447050	1700 cm ²	5000 mL	122 x 499 mm	Seal Cap	Blue	TC Treated	1/12
447060	1700 cm ²	5000 mL	122 x 499 mm	Vent Cap	Blue	TC Treated	1/12

*od=outer diameter (bottom); h=height with closure



Sterile Roller Bottles, Expanded Surface, Ribbed wall, PS

Made from optically clear, virgin, [Polystyrene](#), these bottles have [expanded pleated surface](#) of [2125 cm²](#) and [4250 cm²](#) that increases the growth area by [2.5X](#) and [5X](#) respectively. Ribbed design provides twice the surface area with the same exterior dimensions. These are available either with [seal screw caps](#) or [vent screw caps](#) having 0.2µ PTFE membrane closures that ensure contamination-free cultivation.

Code	Surface Area	Nominal Capacity	Dimensions (od x h)	Cap Style	Cap Colour	Surface	Pack/Case
447121	2125 cm ²	2000 mL	122 x 272 mm	Seal Cap	Blue	TC Treated	1/12
447131	2125 cm ²	2000 mL	122 x 272 mm	Vent Cap	Blue	TC Treated	1/12
447150	4250 cm ²	5000 mL	122 x 499 mm	Seal Cap	Blue	TC Treated	1/12
447160	4250 cm ²	5000 mL	122 x 499 mm	Vent Cap	Blue	TC Treated	1/12

*od=outer diameter (bottom); h=height with closure

Sterile Roller Bottles, Smooth Wall, PET

Made from very stable and excellent gas barrier material [Polyethylene Terephthalate \(PET\)](#). They have a [Smooth surface](#) with a standard growth area of [850 cm²](#). These come in 2 L capacity. Safety screw cap guarantees contamination-free cultivation.

Code	Surface Area	Nominal Capacity	Dimensions (od x h)	Cap Style	Cap Colour	Pack/Case
447041	850 cm ²	2000 mL	122 x 272 mm	Seal Cap	Blue	1/12

*od=outer diameter (bottom); h=height with closure



TARSONS PRODUCTS LIMITED

902, Martin Burn Business Park,
BP-3, Saltlake, Sector - V, Kolkata - 700091
Tel: +91 033 - 35220300, e-mail: info@tarsons.com
Website: www.tarsons.com



TarsonsProductsPvtLtd



tarsons-products-ltd



tarsons.products



TarsonsProduct

